




To begin using HAL you need to make sure module is enable for the system and the appropriate users.

Turn on the HAL option for the system by clicking on the “Hydration Alerts” module on the list in Site Info



Also on the Site Info screen; Opt(1) tab; you need to set your “Weight deficit percentages”

Turn on the HAL option for the applicable users details screen by setting the desired security level for “Hydration Alerts” entry on the module tab...



Note: Settings for the Hydration Alerts (HAL) are applied to the ATS Core, ATS Core Portal, ATS Quick Med Check and ATS Smartphone.